


Science  made smarter


OtoAccess[®] Worklist HL7

EHR and HL7 integration
made easy




Interacoustics

Audiometry


Tympanometry


ABR


OAE


Hearing Aid Fitting


Balance


A clinician seeing
12 patients daily
will typically
free up 1-2 hours
every day with
OtoAccess®
Worklist HL7

Integration
with electronic
healthcare
record (EHR)
systems

Worklist HL7 Integration

Data flow **made** smooth

Database solutions such as OtoAccess® Database provide an easy and secure way to store healthcare data.

However, have you considered how your database fits into your entire workflow? In other words, how do you get data from OtoAccess® Database to your clinic or hospital's electronic health record (EHR) system?

This is where OtoAccess® Worklist HL7 comes into the picture.

What is OtoAccess® Worklist HL7?

This solution automates the communication between EHR systems that support HL7 version 2.x and OtoAccess® Database. What does this mean?

Life without OtoAccess® Worklist HL7

Without OtoAccess® Worklist HL7, your workflow might look like this:

- 1) The patient and appointment are created in your EHR system.
- 2) You look up your patient in your EHR system when he or she arrives.
- 3) You manually create your patient in OtoAccess® Database.
- 4) You perform a certain test, which is stored in OtoAccess® Database.
- 5) You then have two options: either you export a PDF of the test results or you print them.
- 6) At last, you have almost climbed the Mt. Everest of data handling. You either upload the PDF or scan the printout to your EHR system.

This is a long and tedious process! We have spoken with many practitioners, and the process of steps 3), 5) and 6) typically take five to ten minutes. Your time is much better spent doing what you are best at: taking care of your patients.

Life with OtoAccess® Worklist HL7


Let us look at what your workflow would look like with OtoAccess® Worklist HL7:

- 1) Again, your patient and appointment are created in your EHR system.
- 2) Immediately, you can access your patient in OtoAccess® Worklist HL7 and perform your tests.
- 3) Once you are finished, click send. And you are done! Easy, right?

Your job just became a lot easier! No patient creation, no exports and no printing or scanning. This makes your data handling more robust and saves you a lot of time. Think about it. Five to ten minutes saved for every patient. That is a lot of time freed up over the course of a day!

Workflow without OtoAccess® Worklist HL7


Clinicians typically estimate that data-administrative tasks take at least five to ten minutes per patient.



Workflow with OtoAccess® Worklist HL7

The data-administrative tasks are reduced to a few seconds per patient.

Save
5-10 minutes
per patient



Science made smarter

Interacoustics is more than state-of-the-art solutions.

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

Interacoustics.com


Interacoustics A/S

Audiometer Allé 1
5500 Middelfart
Denmark

+45 6371 3555
info@interacoustics.com

interacoustics.com

Full product line support



OtoAccess® Worklist HL7 supports an order-based HL7 v2.x setup in conjunction with OtoAccess® Database. HL7 message types can be configured according to your needs. For details

on configuration, please refer to the OtoAccess® Worklist HL7 compatibility checklist, which will determine if the solution is right for you.

